

## PATIENT POSITIONING DEVICE

### PRODUCT FEATURES

- 1. The Safe, Versatile, Multi-Positioning Device** can position patients into almost any position required by the surgeon and can be stored unobtrusively under the operating table when not in use.
- 2. Eliminate cross infection** by the use of SINGLE USE gel pads fitted to the end of the SurgiPOD Device.
- 3. Faster setup times.** The SurgiPOD Device remains stowed under the operating table between surgical procedures.
- 4. Maintaining Sterility.** The SurgiPOD Device can be repositioned through sterile drapes without compromising sterility.
- 5. Patient Safety.** The SurgiPOD Device has a positive locking tooth assembly in all joints and does not rely on friction between mating surfaces.
- 6. Back Safety.** Theatre staff have no need to lift anaesthetised patients on operating tables, even in emergency repositioning situations.
- 7. More flexibility.** Patients often need extra support or table extensions. This Device provides both.
- 8. Reliability.** Once locked into position this Device will not fail under heavy loads.
- 9. Emergency Repositioning.** In the event of cardiac arrest the patient can be repositioned in seconds, not minutes.

The SurgiPOD Patient Positioning Device is a versatile, multi-position, quick release system that can be released in seconds, not minutes. It allows theatre staff access to the patient at all times and will not obstruct the surgeon from the patient in any way during surgery.



Doc No. A90194 – Sht2 –Page 1-Rev H