

Positioning patients has never been so easy.

“The main benefit is that the patient does not need to be lifted into a particular position to suit the device, because the device has multiple universal joints and hence can be put into any position the Surgeon deems necessary.”

A quote from Professor Berni Einoder, Director of Surgery, Launceston General Hospital, Tasmania



Quick-Connect Coupling allows for other attachments to be fitted, including head rests, arm boards and retractors.

There is no other device on the market which allows such versatility and safety when positioning patients for surgery. **The SurgiPOD Positioning Arm is truly a multi-positioning device.**

The SurgiPOD Patient Positioning System has been specifically designed to:

- **Help reduce the incidence of cross infection through the use of single use gel pads.**
- **Give patients optimal pressure relief.** The single use gel pads can be placed under the patients heels post op and return with the patient back to the ward.
- **Reduce the incidence of skin damage.** The polyurethane membrane won't stick to the patient's skin.
- **Positively lock the patient into the desired position.** The SurgiPOD Device has a quick release mechanism that can easily be operated through the sterile drapes. In the event of cardiac arrest the patient can be quickly repositioned from lateral to supine in seconds.
- **Accommodate all patients.** Irrespective of the patient's size the SurgiPOD Device can deliver optimal strength, flexibility and safety.
- **Fold away and store beneath the operating table, reducing set-up times and storage needs.**
- **Act as one complete system.** Often positioning devices come in many parts and theatre staff don't have time to waste looking for them. The SurgiPOD Device remains on the operating table between procedures.

Doc No. A90194 – Sht2 –Page 2-Rev H