

PELVIC BONE HIP POSITIONER

The SurgiPOD Pelvic Bone Positioning Device is a versatile rigid positioning device that can be setup in under one minute.



PRODUCT FEATURES

1. The SurgiPOD Pelvic Bone positioning device is specially designed for hip navigation surgery holding patients in the lateral position for total hip and hip revision procedures.
2. Eliminate patient to patient cross infections by the use of SINGLE USE Screw on/ Screw off gel pads fitted to the support plate.
3. The specially designed gel pad provides a very even pressure distribution between the patient and the pad resulting in good blood flow under the pad.

4. Patient setup times are reduced because the support plate connects directly to the quick connect coupling on the end of the SurgiPOD positioning arm.
5. The device can be quickly moved away from the operating field to give surgeons and theatre staff access to the patient in emergency repositioning situations.
6. The device can be adjusted through the sterile drapes.
7. The device can be autoclaved.



Doc No. A90194 – Sht5 –Page 2-Rev I